MANUALE ISTRUZIONI  
OPERATOR'S HANDBOOK  
MANUEL D'UTILISATION  
GEBRAUCHSANWEISUNG

# IDROCICLONE ISO 6

# IDROCICLONE ISO 6

## INDICE

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# HYDROCYCLONE ISO 6

## INDEX

Introduction.....	12
Signals and warnings .....	12
Field of use .....	12
Working .....	12
Technical features .....	13
Connection .....	14
Directions for use .....	14
Installation .....	16
Mounting .....	16
Notices.....	17
General warnings and biological risk .....	17
Routine maintenance entrusted to specially instructed surgery staff .....	17
Faults report and possible solutions .....	18
Maintenance and periodic checks reserved for authorized technicians .....	18
Warranty .....	19
Transport of second-hand appliances .....	19
Important notices .....	19
Drawings.....	40

- *Introduction*
- *Signals and warnings*
- *Field of use*
- *Working*

- **Introduction**

The following presentation is intended to illustrate the equipment and the appliance in order to clarify their working and maintenance. It also informs users and technicians of possible dangers and about the precautions to be taken for accident prevention.

- **Signals and warnings**

- *Electrical shock risk : also 230 V  $\sim$  can be lethal.*

- *Biological danger, risk of infections from epidemic diseases.*

- *General danger sign.*

- *High temperature.*

- *Compulsory direction of flow or of rotation.*

Signs cannot fully express danger warnings and compulsory notices, therefore it is necessary that the user reads the warnings and keeps them in due consideration. Failure to observe a danger sign or warning may harm operators or damage the machine.

Safety devices must not be removed. The structure of the appliance/equipment or their functioning must never be modified.

- **Field of use (only for the German market)**

According to the annex no.50 of the General Framework of the Administrative Standard regarding the minimum requirements for conveying waste water coming from dental surgeries to the drainage system, it is necessary to install amalgam separators to treat such water containing amalgam.

- **Working (draw. 1 page 40)**

Inside the centrifugal separator air is separated from aspirated liquid, then after passing through the inlet passage (1) the aspirated liquid reaches the evenly spaced holes (2) which spread the fluid uniformly inside the cone and give it a rotary descending motion.

When the fluid reaches the bottom of the cone (4), inversion of the flow rotation and direction takes place; at that moment of zero speed amalgam particles fall into the bottom of the container (5), whereas the liquid rises out of the middle of the cone.

The re-circulation valves 6 and 7 open and close the passage for the drainage of the amalgam depurated liquids.

During the phase A the system has not reached the fixed flow, therefore the liquid is conveyed back to the Hydrocyclone; in the phase B, with the suitable flow, separation is obtained and the processed liquid can be conveyed to the waste (8).

- *Technical features*

Name	HYDROCYCLONE ISO 6
Classification according to ISO 11143	Type 4
Maximum flow according to ISO 11143	5,5 l/min.
Degree of separation according to ISO 11143 with 5,5 l/ min. flow	> 98%
Degree of separation according to ISO 11143 with 0,5 l/ min. flow	> 98%
Amalgam container capacity	115 ml
Room temperature	max. + 40 °C, min. + 5 °C
Control panel voltage drive	12 V DC
Amalgam level sensor voltage drive	12 V DC
Dimensions	H= 315 mm , W= 205 mm , D= 165 mm

Hydrocyclone ISO 6 is supplied with the “Separatore Centrifugo Micro-Smart” (Centrifugal Separator Micro-Smart) of which you will find some details listed below.

Name	“SEPARATORE CENTRIFUGO MICRO-SMART”
Electric power supply	230 – 400 V 3 ~ 75 Hz 3-1,7 A
Maximum allowed negative pressure	300 mbar

- *Connection*
- *Directions for use*

- **Connection**

- *Electrical shock risk : also 230 V ~ can be lethal.*

- *General danger sign.*

- *Biological danger, risk of infections from epidemic diseases.*

- *Compulsory direction of flow or of rotation.*

This section is intended to illustrate the connections between the Hydrocyclone ISO 6 and the “Separatore Centrifugo Micro-Smart” (Centrifugal Separator Micro-Smart).

- **Hydraulic connection (draw. 3 page 42)**

All the connections must be properly fixed and secured using a mechanical securing system like a pipe clamp.

For the connections use only original hoses. Replace them every two years at the latest. Do not use hoses with a diameter different from the one of the hoses originally supplied with the Hydrocyclone.

- |                                                  |                           |
|--------------------------------------------------|---------------------------|
| • Hose inlet connection (11) to the Hydrocyclone | Hose diameter d = 11,7 mm |
| • Re-circulation hose connection (12)            | Hose diameter d = 11,7 mm |
| • Connection for drainage to the waste (8)       | Hose diameter d = 20 mm   |

- **Electric connection (draw. 4 page 43)**

- |                                               |                      |
|-----------------------------------------------|----------------------|
| • Control keypad cabling                      | Cable with RJ 9 plug |
| • Amalgam level sensor (13) cabling on AC 86C | Cable with RJ 9 plug |

- **Directions for use**

**Description of the electronic amalgam level control and alarm system**

Hydrocyclone ISO 6 can be controlled using the electronic control panel AC 100.

With this contact it is possible to place a signal remotely.

**Description of the electronic amalgam level control and alarm system (draw. 6 page 45)**

Through a pair of photodiodes the amalgam level sensor (13) fitted aside the amalgam collection container (5) detects the amalgam level inside the container: the text “Report Cod. I14 ” on the display (15) warns that the amalgam level has reached 95%; press the button ESC (16) to cancel the alarm from the display and to visualize other parameters, remember that the control panel keeps the alarm in its record and that the alarm will be displayed again upon switching on the machine.

The Hydrocyclone keeps on running regularly. When the amalgam level reaches 100%, the text “Report Cod. I15” appears on the display (15). The Hydrocyclone keeps on running until the machine is switched off. It won't be possible to start the machine again unless the amalgam collection container (5) is replaced.

Every time one of the above mentioned alarms occurs, a volt-free (clean) contact closes on the circuit board (terminals 1-2) (draw. 4 page 43). With this contact it is possible to place a signal remotely.

The contact remains closed and the alarm remains active until the amalgam collection container is replaced and the machine is switched off.

- **Replacement of the amalgam collection container (draw. 7 page 46)**

- **Electrical shock risk : also 230 V  $\sim$  can be lethal.**

- **General danger sign.**

- **Biological danger, risk of infections from epidemic diseases.**

- **Compulsory direction of flow or of rotation.**

Before intervening, disconnect from the mains the machine on which the Hydrocyclone is mounted, wear disposable gloves, mask and goggles to protect from possible splashes.

Slowly unscrew the handle (17) as shown in the drawing. Support the container while lowering it. Continue to turn the handle until it reaches its limit-stop. Remove the container keeping it horizontally and fill it with a suitable disinfectant (refer to the importer for suggestions regarding the most suitable product). Open the valve (18) placed on the cover (19), treat the “O” ring with a thin layer of Lubri-Jet, position the cover on the container and close the valve (18).

The container is made of unbreakable material and is fitted with a hermetically sealed lid, since the contained material is toxic and contaminating for epidemic diseases (we recommend to handle it with care).

To receive a list of authorized companies with licence for transport and disposal of this kind of waste and for further information regarding the disposal, contact trading associations or local competent authorities.

Position a new container in the seat of the original one. Screw the handle (17) until the rim of the container comes into contact with the cone, which should not require force.

## • Installation

## • Mounting

### • Installation (draw. 2 page 41)

- *Electrical shock risk : also 230 V  $\sim$  can be lethal.*

- *General danger sign.*

- *Biological danger, risk of infections from epidemic diseases.*

- *Compulsory direction of flow or of rotation.*

Install the Hydrocyclone below the separator in the vertical position as follows:

Mount the bracket (10) without disassembling the Hydrocyclone, as shown in the drawing. The Hydrocyclone must be easily accessible from every side. Leave a free space of at least 10 cm under the amalgam container (5) so it can be easily removed.

The filter (14) fitted to the centrifugal separator has a particle size not larger than 1.2 mm, protecting the Hydrocyclone from larger particles which may damage it.

### • Mounting (draw. 5 page 44)

- *Electrical shock risk : also 230 V  $\sim$  can be lethal.*

- *General danger sign.*

- *Biological danger, risk of infections from epidemic diseases.*

- *Compulsory direction of flow or of rotation.*

#### **Distances**

The maximum distance does not depend on the Hydrocyclone but on the suction pump, we therefore recommend that you refer to the instructions of the suction pump supplied by the manufacturer.

#### **Positioning**

The Hydrocyclone can be installed in the following layouts:

- 1) installation on the same level of the dental units
- 2) installation on the lower level

The Hydrocyclone must be easily accessible from all sides. Leave a free space on the top and on the front of the unit to facilitate maintenance procedures.

Liquids drain due to fall, therefore we recommend to provide a fall towards the sewage to facilitate the draining. It is possible to place a small trap, provided that it does not obstruct the draining.

- *Notices*
- *General warnings and biological risk*
- *Routine maintenance entrusted to specially instructed surgery staff*

- **Notices**

After installation, the installer should train the surgery staff with practical demonstrations for amalgam container replacement and filter cleaning on new and therefore not contaminated machines.

Installation, connection, mounting, maintenance and other work on the amalgam separator must be entrusted to technicians who have been specially trained in the dental field and who are authorized by the manufacturer. The aspirated liquid is contaminating and infected, therefore we recommend to take all the necessary precautions to prevent personal and environmental contamination. Bad running can also lead to contamination.

- **General warnings and biological risk**

- Before starting any servicing operation on appliances which have been used, carry out a series of rinses with a suitable disinfectant.

- Disconnect the centrifugal separator from the mains supply and padlock the isolating switch, if present.

- Before starting any maintenance work, wear gloves, goggles, mask and disposable overalls.

The centrifugal separator pressurises aspirated liquid and forces it inside the Hydrocyclone. In case of a failure of a pipe or of a component under pressure biological risk is more serious.

- **Routine maintenance entrusted to specially instructed surgery staff (draw. 3 page 42)**

- *Electrical shock risk : also 230 V  $\sim$  can be lethal.*

- *General danger sign.*

- *Biological danger, risk of infections from epidemic diseases.*

- *Compulsory direction of flow or of rotation.*

1) After each surgical operation or particularly long operations, aspirate abundant clean water to rinse the suction piping.

2) At the end of each working day, clean and disinfect the appliance using a specific product for dental aspirators and amalgam separators. We suggest, for example, Puli-Jet plus (for dilution with water, refer to our directions). Do not use aggressive or foaming products.

3) The dirty filter (14) must be replaced every day. Handle it with care to prevent contamination and the return of solid trapped particles into the filter. Any material trapped in the filter must be collected and disposed of according to the regulations in force.

4) We recommend placing one disinfectant antifoam tablet in the dental unit secretion filter every morning before starting the aspiration.

- *Faults report and possible solutions*
- *Maintenance and periodic checks reserved for authorized technicians*

5) In case the appliance is treated with products which are different from the ones recommended by the importer, the correct working of the equipment might be compromised and faults could occur.

6) Check that there are no liquid leaks. In case of working anomalies or alarms, contact the technician.

- ***Faults report and possible solutions (draw. 6 page 45)***

The amalgam container filling levels are the only alarms directly linked to the Hydrocyclone.

Below is information useful to inform the reader about possible problems that may cause the aspiration to stop and that may compromise the working of the “Separatore Centrifugo Micro-Smart” (Centrifugal Separator Micro-Smart) for the Hydrocyclone.

1) If the display (15) is off it means that there is no power supply: check that the machine switch and the master switch are on, in case they are on but the machine is still off, contact the technician.

2) The text “ Report Cod. S09, S10 or S11 “ indicates that the AC 100 amperometric control system has detected an anomaly in the centrifugal separator current absorption.

Such an anomaly could have different causes: a mechanical blockage or a drain obstruction (8), (draw. 3 page 42) of the Hydrocyclone or a temporary voltage decrease.

Disconnect the appliance from the mains for a few seconds to reset the alarms and to start the appliance again. In case the alarm appears again, contact the technician.

- ***Maintenance and periodic checks reserved for authorized technicians***

- ***Electrical shock risk : also 230 V ~ can be lethal.***

- ***General danger sign.***

- ***Biological danger, risk of infections from epidemic diseases.***

- ***Compulsory direction of flow or of rotation.***

**Yearly maintenance (draw. 8 page 47)**

1) Clean the suction piping by aspirating at least 2 l. of a solution prepared with a detergent-disinfectant for dental aspirators.

2) Remove the amalgam collection container (5) as per the instructions of the paragraph “Replacement of the amalgam collection container” and check the working of the level sensor (13) inserting a metallic piece between the two photodiodes. The Hydrocyclone control panel or the display will show a visual alarm to indicate that the amalgam level has reached 95% or 100%.

3) Remove the Hydrocyclone (20) upper part, clean all its elements, especially the five holes (2), check the conditions of the ‘O’ rings and lubricate them with Lubri-Jet before mounting them again.

- *Warranty*
- *Transport of second-hand appliances*
- *Important notices*

- 4) Check and clean, if necessary, the inside of the cone, in particular the draining hole.
- 5) Check that the connection piping is not leaking and that it is fixed. If in doubt, replace it. The connection piping must be replaced every two years.
- 6) Check carefully that there are no leaks.

#### **Five-yearly maintenance (draw. 9 page 48)**

\*According to the annex no.50 of the General Framework of the Administrative Standard regarding waste water, amalgam separation systems must be checked every five years at the latest to make sure that they still comply with the standard.

After yearly maintenance checks have been carried out it is necessary

- 1) to remove and clean the valves (6) and (7) and all the elements of the re-circulation system, to check the conditions of the internal elements and replace them, to check and replace, if necessary, all the 'O' rings lubricating them with Lubri-Jet before mounting them again;
- 2) to make sure that connections, installation and mounting comply with the instructions of this manual.

#### • **Warranty**

The appliance is guaranteed for one year from the date of sale, provided that the warranty slip is returned to the manufacturer with date of sale and customer's name. Warranty and manufacturer's liability cease in case the appliance and the equipment are tampered with operations of any kind carried out by people who are not authorized by the manufacturer.

#### • **Transport of second-hand appliances**

- **Biological danger, risk of infections from epidemic diseases.**

Before packing, clean and disinfect with a suitable disinfectant. Remove the liquids and residual debris which are inside the appliance, place it in a polyethylene bag, seal and pack it in packaging constructed of 3-layer corrugated board.

#### • **Important notices**

- The manufacturer, concessionaries, agents and authorized technicians are at customers' disposal for advice and assistance, to supply literature, spare parts and anything useful.
- The manufacturer reserves the right to modify the products for improvements, for technical, normative and functional reasons or for problems due to the availability of products or semi-finished products, without prior notice.
- Our updated manuals are available on the web site [www.cattani.it](http://www.cattani.it). We recommend they are consulted especially for updates about safety.

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