

asametro

CATTANI S.p.A. 43100 Parma (Italy) via G. Natta 6/a - tel. 0521-607604 - fax 0521-607628
ESAM S.p.A. 43100 Parma (Italy) via G. Natta 4/a - tel. 0521-607613 - fax 0521-607628

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ASAMETER

The asameter is a variable space appliance based on the principle of Bernoulli.

Essentially the appliance is made of a polysulfone tube with afloat; the tube is conic inside and gauged with very close tolerances.

When a fluid passes through the appliance the afloat moves inside the tube and stops at the level corresponding to the flow rate.

The values of the flow rates can be read on a scale on the tube.

The value is given by the scale at the point reached by the upper part of the afloat.

HOW TO TEST OUR ASPIRATORS AND CENTRALIZED PLANTS WITH OUR ASAMETER «P13».

How to test our aspirators and centralized plants with our Asameter P13.

The Asameter must be in the upright position and the passage of air should be free of any obstacle. The «A» reduction has been designed to check the aspiration on the small terminal (item 1115) of our tip support; the «B» reduction is to be inserted into the large terminal (item 1114); the tube-holder «C» is useful to check the aspiration directly on the 30 Ø pipe.

In case the flow test is to be done directly on the aspiration unit, the suitable Asameter, with the proper capacity must be requested, making reference to the working diagrams on the various technical sheets.

Please note that it is impossible to test simultaneously the large and the small terminal on the same tip support: the two hoses start very closely from the same manifold and the flow passes through the large one rather than the small. We are sure that 250/300 litres of air per minute through the large terminal and 120/150 through the small one are enough to service properly any operation; higher flow rates could disturb patients suffering from sensitivity at the teeth collars or on other places not anaesthetized.

P.S. The tolerances in our Asameter is about 4%. The results of the Asameter could be modified by the following factors:

- 1) Change of pressure on the measurement point
- 2) Change of temperature at the measurement point.

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