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## **INDEX**

**Page**

— REFERENCES .....	11
— INTRODUCTION AND FUNCTIONING HYDROSEPARATOR .....	12
— INTRODUCTION AND FUNCTIONING HYDROCYCLONE .....	13
— HOW TO REPLACE THE AMALGAM CONTAINER .....	13
— MAINTENANCE AND CLEANING .....	14
— EXTRAORDINARY EVENTS .....	15
— NOTICE .....	15
— WARRANTY .....	15
— TRANSPORT OF SECOND-HAND APPLIANCES .....	15
— HYDROSEPARATOR TYPE DIN BASIC DRAWING .....	41
— HYDROSEPARATOR TYPE ISO BASIC DRAWING .....	42
— HYDROCYCLONE TYPE DIN BASIC DRAWING .....	43
— HYDROCYCLONE TYPE ISO BASIC DRAWING .....	44
— HYDROCYCLONE TYPE ISO BASIC DRAWING FOR TURBO-JET MODULAR .....	45
— SPITTOON DRAINING BASIC DRAWING .....	46

## REFERENCES

Tables page 41 42 43 44 45 46

### HYDROSEPARATOR

- 1 - *Mini-Separator*
- 2 - *Draining unit*
- 3 - *Hydrocyclone*
- 4 - *Mini-Separator debris inlet*
- 5 - *Overflow probes*
- 6 - *Air outlet*
- 7 - *Hydrocyclone impeller*
- 8 - *Hydrocyclone debris inlet*
- 9 - *Recycling electrovalve*
- 10 - *Draining electrovalve*
- 11 - *Hydrocyclone draining electrovalve*
- 12 - *Hydrocyclone draining switch*
- 13 - *Hydrocyclone reset switch*
- 14 - *Beeper reset*
- 15 - *Power-on signal*
- 16 - *Amalgam level 95%*
- 17 - *Amalgam level 100%*
- 18 - *Debris maximum level*
- 19 - *Spittoon electrovalve*
- 20 - *Spittoon draining unit*
- 21 - *Amalgam container*
- 22 - *Hydrocyclone small motor*
- 23 - *Amalgam container handle*
- 24 - *Amalgam container cover*
- 25 - *Amalgam container valve*
- 26 - *Electronic P.B.*
- 27 - *Electropneumatic valve*
- 28 - *Hydroseparator box*
- 29 - *Box cover*
- 30 - *Fixing bracket*
- 31 - *Control panel*
- 32 - *Spittoon draining unit fitting*
- 33 - *Liquid drainage tap*
- 34 - *Switch for emptying the Mini-Separator*

### HYDROCYCLONE

- Collecting tank*
- Draining unit*
- Hydrocyclone*
- Collecting tank liquid inlet*
- Overflow probes*
- Air outlet*
- Hydrocyclone impeller*
- Hydrocyclone liquid inlet*
- Recycling electrovalve*
- Draining electrovalve*
- Hydrocyclone draining electrovalve*
- Hydrocyclone draining switch*
- Hydrocyclone reset switch*
- Beeper reset*
- Power-on signal*
- Amalgam level 95%*
- Amalgam level 100%*
- Debris maximum level*
- Spittoon electrovalve*
- Spittoon draining unit*
- Amalgam container*
- Hydrocyclone small motor*
- Amalgam container handle*
- Amalgam container cover*
- Amalgam container valve*
- Electric circuit*
- Spittoon draining unit fitting*
- Liquid drainage tap*

# **HYDROSEPARATOR** (for air systems)

## **DIN and ISO type**

(all listed general features, not mentioned for a single model, refer to both DIN and ISO type).

## **Introduction** (part. page 41 - 42 - 46)

The amalgam separator for air aspiration systems - **Hydroseparator** - has been obtained by combining the **Hydrocyclone (3)** to the **Mini-Separator (1)**.

The Hydroseparator is produced in four types: without box - to be installed inside the dental unit, for single surgery and for centralized systems and with box - to be installed outside but next to the dental unit for single surgery and for centralized systems. Two other models are available: for semi-wet and wet systems; in such systems the amalgam separator must be installed in the plant room without the Mini-Separator and it is then called Hydrocyclone. **The DIN-type Hydrocyclone has been approved by TÜV for a 97% separation rate; the ISO type has been approved by TÜV for a 96,6% separation rate.** The Hydroseparator works always in presence of liquids and is full at all times even when the system stops (see more details in the "functioning" paragraph). The continuous presence of liquids and a proper and methodical maintenance, carried out with the suggested materials, help keep the appliance clean and in a good working conditions. **The functioning of the Hydroseparator and Hydrocyclone is totally automatic. In the DIN version functioning is controlled by the low tension electronic circuit (26) which in its turn is operated by the electrical circuit of the aspiration system. Functioning is automatic for the ISO version as well; upon request it can be equipped with electronic sensor.**

## **Functioning** (det. page 41 - 42)

When one of the hoses is lifted from its seat, the Electropneumatic valve (27)\*\* opens, the aspiration unit is switched on and aspiration starts. **Aspirated fluids enter the Mini-Separator (1):** air gathers in the middle and is drawn in by aspiration power, liquids to the walls and drop to the bottom. When level of liquids reaches the second probe of the level control unit (5), the draining unit (2) is switched on, the valve opens and **liquids are conveyed into the Hydrocyclone.**

The Hydrocyclone impeller (7)\* **starts spinning making liquids and amalgam rotate inside the cyclone: when liquids reach the bottom of the cyclone, rotating direction and flow are reversed. For one second the liquids are still and the amalgam particles drop down into the amalgam container (21);** then the liquids flow back to the top (from the centre of the cone) and are drained out. Liquids in the spittoon go through the draining unit (20) and consequently aspiration starts and the Hydroseparator with it. Spittoon waste is also purified from amalgam. In order to avoid overworking the Hydroseparator, spittoon rinsing must be timed at max. 20 seconds.

\* Only for DIN type

\*\* Only in centralized plants

# **HYDROCYCLONE** (for semi-wet and systems)

## **DIN and ISO type**

### **Introduction** (det. page 43 - 44 - 45 - 46)

**The Hydrocyclone is an amalgam separator for semi-wet and wet systems.** Accessories for the separator are different for air, wet and semi-wet systems. For wet and semi-wet systems the Hydrocyclone and its accessories are to be located in the plant room, while **the spittoon draining unit (20) with its electrovalve (19) will be installed inside the dental unit** next to the spittoon.

The spittoon outlet is to be connected to the upper part of the draining unit (20A), outlet (20B) is to be connected to the aspiration piping coming from the tip support through adapter (32). Moreover, **the spittoon unit (20) shall be at the surgery staff's reach to allow for replacement of the filter** inside (20C).

When installed together with big Separator-Tanks, the Hydrocyclone (ISO or DIN) must be located between the draining pump and the waste outlet. All Turbo-Jet Modular units can be combined with an ISO amalgam separator and can be purchased with or without separator. Wet-ring pumps Pal 22 and 50 can be combined with a DIN amalgam separator with reservoir, installed between the pumps and the waste outlet, or with a centrifugal separator and ISO amalgam separator installed between the chair and Pal 14-22-38-50.

The DIN Hydrocyclone draining flow is 3.5 l/min; the ISO type draining flow is 5.5 l/min. During washing and disinfection of the plant there is a higher risk of interruptions owing to overflow; however interruptions are momentary and limited to washing time. The collecting tank (1) is suitable for installation of two Hydrocyclones, allowing for a 7 l/min flow. In semi-wet plants with big Separator-Tanks and Turbo-Jet modular you need to time spittoons at 20 seconds and wait a few seconds in between washing sessions: this will enable the draining pump and centrifugal Separator to drain aspirated liquid. With the "smart bucket" washing can be performed without interruptions.

### **Functioning** (det. page 43 - 44 - 45)

Some filters for larger amalgam particles are installed before the amalgam separator. In the separator unit either the draining pump or the centrifugal separator drive pre-filtered aspirated liquids into the cyclone (ISO or DIN); separation is performed inside the cone of the cyclone. For further details on the functioning of the ISO and DIN Hydrocyclone see paragraph "FUNCTIONING" of the Hydroseparator as the items are the same.

### **How to replace the amalgam container** (det. page 41 - 42 - 43 - 44 - 45)

**On the side of the amalgam container (21) two photodiodes sense the level of amalgam** collected inside the container: a yellow warning light (16) and a beeper indicate that the container is **95%** full. Press switch (14)\* and the beeper will stop, while the Hydroseparator keeps functioning normally and the warning light (16) is still on. Indicator (17) and the beeper will then warn that the level of amalgam has reached **100%** of the container capacity. By pressing switch (14) you reset the beeper, while the Hydroseparator keeps functioning until the system is turned off. The system will not work until the container has been replaced. The ISO type allows for

visual checking of the container. **Visual checking is to be carried out very often, every 30 days max. The amalgam level must not exceed the mark. Prior to removing the container, wear disposable gloves, mask and goggles in order to protect yourself against accidental splashing.**

Keep switch (12)\* pressed for 5 seconds in order to drain residual liquid from the Hydrocyclone and the amalgam container. For ISO version, drainage of residual liquid is obtained by opening tap (33).

Hold the container with your hand to prevent it from falling, release the handle lifting it upwards, **remove the container, and fill it up with pure Puli-Jet; open valve (25) on the cover (24) of the container, place cover on the container and close valve.** The container is made of shock-resistant material and is sealed; the contents are **toxic, harmful and contaminating** for epidemics and shall be handled with care.

The waste amalgam container can also be sterilized in autoclave with the liquid cycle: fill with Solid-Steril only prior to sterilization. Therefore you can label the container as "Container and contents are sterilized in autoclave".

In order to identify the authorized institutions for disposal of waste in your area, contact your local authorities. Further information can be obtained from your national dentists' associations. Abiding by the regulations in force, **the container shall be delivered to a firm or institute which is authorized to transport, dispose or recycle such substances. A new container must be fitted to the system** and the handle properly tightened onto the container. In order to start the system over again keep switch (13)\* pressed for 5 seconds and the separator will start functioning. For ISO version just close tap (33).

When the container has filled up to 95% of its capacity, you can wait for maximum level; however you run the risk of a stop, as in case you should turn off the appliance you could not turn it back on. **Once reached 95%, we suggest you to clean the appliance with Puli-Jet prior to replacing the container.** This way you will considerably limit the risks of contamination and forced stop.

## **Maintenance and cleaning (det. page 41 - 42 - 46)**

Liquids coming from the spittoon and the tip support are to be **filtered through a net filter capable of retaining any particle greater than 0.8 mm.**

The spittoon draining unit (20) is provided with another filter (20C) that needs to be checked regularly. Such filter is subsidiary and cannot be used instead of the spittoon filter which does not always meet the safety requirements of the Hydroseparator (retaining particles greater than 8 mm).

With the Hydroseparator and the Hydrocyclone we suggest the use of sanitizing antifoaming agents (see instructions in the package); the antifoaming agent inhibits foam build up and prevents a possible accidental stop of the appliance, meanwhile reducing bacterial rate inside the Hydroseparator and piping.

**At the end of each working day, it is important to aspirate a Puli-Jet solution in hot water (50 °C) in all types of systems, air, semi-wet and wet systems.** The solution must be prepared following the instructions on the bottle and proceed as follows: fit tip No 11, 20 or 17 to the large terminal and No 10, 17 or 21 to the small terminal; dip one tip at a time into the solution (do not dip the terminal) and aspirate until suction power is reduced owing to the liquid, remove the tip from the solution and lift it up so that the liquid moves quickly from the hose to the manifold and on to the Mini-Separator (or to the piping in case of wet or semi-wet systems). Repeat 2 or 3 times for each terminal.

\* **Only for DIN type and ISO type with electronic level control (ISO type does not have switch 12 and 13; in order to empty the Mini-Separator press switch 34)**

By fitting the terminals to the plungers of Cattani's 'smart bucket' the above operation is carried out automatically. **Regular use of Puli-Jet assures correct functioning and preservation of the Hydroseparator, while that of other products could jeopardize it.** The manufacturer disclaims all responsibility and shall not be liable for goods under warranty if they have been treated with other products.

## **Extraordinary events\* (det. page 41)**

In case of a blackout while the Hydroseparator is in service and the level control probes are dipped in liquid, when the power supply is restored the beeper will go off, the warning light (18) will indicate maximum level: press switch (14) to reset the appliance.

The same will take place in case of a power failure while replacing the amalgam container: in this case, prior to resetting the appliance with switch (14), make sure you have fit a new amalgam container.

## **Notice**

The retailer or installer will instruct the dental office staff with practical tests on new uncontaminated appliances. In case the Hydroseparator should stop owing to unknown failures or reasons, do not work on the appliance. **All repairs and part replacements are meant to be carried out by authorized technicians who are provided with original spare parts.** Installation of the amalgam separator is to be carried out by technicians of the dental field, authorized by the manufacturer.

All three versions of amalgam separator - air, wet and semi-wet versions - are supplied pre-assembled and cabled.

**Aspirated liquid is always contaminated and infectious, therefore great care is to be taken to protect users and the environment from contamination. Mal-functioning may also produce contamination, therefore be sure you address only to skilled qualified technicians who have been trained and conferred a certificate by the manufacturer.**

Any modification to the appliance shall be approved by the manufacturer. **For any situation not listed in these instructions, the manufacturer, distributors and agents are willing to provide the necessary information.**

## **Warranty**

- The appliances are guaranteed for one year from date of sale, provided that guarantee card addressed to manufacturer is returned to the manufacturer reporting date of sale, retailer's stamp and customer's name.
- Guarantee and manufacturer liability cease in case the appliances and plants are found tampered by any kind of action performed by unauthorised people.

\* Only for DIN type

## ***Transport of second-hand appliances***

*Prior to packing, clean and sanitize with Puli-Jet (see "Cleaning and Maintenance" section).*

*Place unit into a polyethylene bag, seal and pack in 3-layer corrugated board.*

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**CATTANI S.p.A.**
**ESAM S.p.A.**